

GOT (AST, SGOT)

IFCC

(Glutamate Oxaloacetate Transaminase)

Diagnostic reagent for determination of AST activity.

Liquid. Dual Reagents. Store at 2°C - 8°C. For in Vitro Diagnostic Use. Do not freeze

Ref No	Pack						
A2210	5 x 100 mL	DM2210	960 Tests	BY2210	7670 Tests	K2211	4000 Tests
A2211	5 x 50 mL	R2210	2885 Tests	BY2211	5114 Tests	M2210	2512 Tests
A2212	5 x 25 mL	R2211	962 Tests	N2210	800 Tests	M2211	1675 Tests
T2210	4000 Tests	S2210	3077 Tests	N2211	400 Tests	L2210	3750 Tests
T2211	2121 Tests	S2211	1477 Tests	K2210	6000 Tests	L2211	2000 Tests

INTENDED USE

The test is applied for the quantitative determination of AST activity in serum and plasma.

TEST PRINCIPLE

The enzyme aspartate aminotransferase (EC 2.6.1.1; L-Aspartate: 2-Oxoglutarate Aminotransferase, AST or AspAT; Glutamate Oxaloacetate Transaminase, GOT) catalyzes the transaminase reaction between L-Aspartate and 2-Oxoglutarate. The 2-Oxalacetate formed, is reduced to malate in the presence of MDH. As the reactions proceed, NADH is oxidized to NAD. The disappearance of NADH per unit time is followed by measuring the decrease in absorbance at 340 nm. The present method has been made according to IFCC (2002).

TEST PARAMETERS

Method: UV, Kinetic, Decreasing Reaction,

mod. IFCC

Wavelength: 340 nm Temperature: 37°C

Sample : Serum, EDTA-Plasma,

heparinized Plasma

Linearity : - 440 U/L (340 nm, 37°C)

REAGENT COMPOSITION

 $\begin{array}{lll} \mbox{Tris buffer} & \leq 90 \mbox{ mM pH 7.65}, \\ \mbox{L-aspartate} & \leq 250 \mbox{ mM}, \\ \mbox{2-Oxoglutarate} & \leq 14 \mbox{ mM}, \\ \mbox{NADH} & \leq 0.18 \mbox{ mM}, \\ \mbox{MDH} & \geq 600 \mbox{ U/L}, \\ \mbox{LDH} & \geq 900 \mbox{U/L}. \\ \end{array}$

REAGENT PREPARATION

Mix 4 parts of Reagent 1 (Buffer, Enzymes) with 1 part of Reagent 2 (Substrate). For example: 4 ml Reagent 1 and 1 ml Reagent 2.

Working reagents are stable at 2-8°C in case of closed vials and avoiding contamination after preparation.

For manual working procedures; if working reagent will be used; shake the Reagent 2 vial gently before pouring its contents into the Reagent 1 bottle. It is advisable to wash the Reagent 2 vial with a small volume of the prepared mixture in order to completely rinse the vial and avoid any losses.

REAGENT STABILITY AND STORAGE

Stability of working reagent is 30 days at 2-8°C, away from light sources. Protect from light. Note expiration date on the label.

Close immediately after use. Avoid contamination of the opened reagents.

Once opened vials are stable minimum 30 days at 2-8°C at optimum conditions. On board stability is strongly related to auto analyzers cooling specification and carry-over values.

SAMPLE

Serum and plasma are collected by standard procedures.

GOT is stable for 4 days at 2-8°C and 1 month at - 20°C.

TEST PROCEDURE

Sample Start

There have many ready application procedures dedicated to different kind of photometers and

Rev: V2.3 Date: 01.13 GOT (AST, GOT) Page 1/3



ready manual working process can be supplied on request.

There have many ready application procedures dedicated to different kind of biochemistry auto analyzers can be supplied on request.

Substrate Start

There have many ready application procedures dedicated to different kind of biochemistry auto analyzers can be supplied on request.

CALCULATION

Perform calculation in units per litre, multiplying the ΔA /min by the factor as it is indicated.

Calculation in U/L: Δ A/min x 1746 Activity in μ kat/I: U/L x 0.0167 = μ kat/L

REFERENCE INTERVAL (NORMAL VALUES) (Based on CLSI C28-P Document)*

Men : < 35 U/L ($< 0.58 \mu kat/L$) Women: < 31 U/L ($< 0.52 \mu kat/L$)

*It is recommended that each laboratory establish its own reference range.

QUALITY CONTROL AND CALIBRATION

All control sera with GOT values determined by this method can be used. We recommend:

"ARCON N", Assayed Control Serum Normal Cat.No. A3910

"ARCON P", Assayed Control Serum Abnormal Cat.No. A3920

The use of a GOT Calibrator (for automated Systems) is optional. We recommend: ARCAL AUTO Calibrator **Cat. No.** A39050

*Calibration Stability: It is strongly depend of application to auto analyzers and auto analyzers specification. Calibration stability is 30 days.

*Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

Quality control is recommended every morning. Calibration is not recommended if QC control values are acceptable. Reagent should be calibrated after lot changes.

PERFORMANCE CHARACTERISTICS

Low linearity:

High Linearity: The method is linear up to 440 U/L.

If a ΔA /min of 0.200 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Linearity may considerably vary depending on the instrument used.

Precision Studies (Based on CLSI EP5 Doc.):

Repeatibility (within run) (Intra-assay)

- CO. C			
Mean conc.	SD	CV	n
38.88 U/L	0.56	1.40%	10
132.50 U/L	0.96	0.70%	10

Reproducibility (run to run) (Inter-assay)

Mean conc.	SD	CV	n
37.87 U/L	0.46	1.20%	20
134.48 U/L	1.48	1.10%	20

Sensitivity (LOD) (Based on CLSI EP17 document): The limit of detection is 0.47 U/L.

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

Interferences: No interference was observed by the presence of:

 $\begin{array}{ll} \mbox{Hemoglobin} & \leq 400 \mbox{ mg/dL} \\ \mbox{Bilirubin} & \leq 25 \mbox{ mg/dL} \\ \mbox{Lipids} & \leq 500 \mbox{ mg/dL} \end{array}$

Methods comparison: A comparison between Archem and a commercially available product gave the following results:

GOT Archem = x GOT competitor = y n = 112 y = 0.986x + 1.636 U/I r2 = 0.997

These performance characteristics have been obtained using an analyzer. Results may vary if a different instrument or a manual procedure is used.

NOTES

- For in vitro diagnostic use only. Do not pipette by mouth. Avoid contact with skin and mucous membranes.
- All the calibrators, controls and some reagents must be considered as human & animal sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.
- Material safety data sheet will be supplied on request.
- Exercise the normal precautions required for handling laboratory reagents.
- After measurements are taken, reagent bottles should cap and kept at 2-8°C. Caps of the reagents bottles cannot be used between two different kind of reagent and between R1 & R2.
- Reagents with different lot numbers should not be interchanged or mixed.
- The linearity limit depends on the sample to reagent ratio.

Rev: V2.3 Date: 01.13 GOT (AST, GOT) Page 2 / 3



PRECAUTIONS AND WASTE DISPOSAL

This product is made to be used in professional laboratories and by professional operators. Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

R36/38: Irritating to eyes and skin.

S20/21: When using, do not eat, drink or smoke.

S26 : In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S28 : After contact with skin wash immediately with plenty of water.

S36/37/39: Wear suitable protective clothing, gloves and eye/face protection.

S45 : In case of accident or if you feel unwell, seek medical advice immediately.

S56 : Dispose of this material and its container at hazardous or special waste collection point.

S57 : Use appropriate container to avoid environmental contamination.

S61 : Avoid release in environment. Refer to special instructions/safety data sheets.

Please consult local regulations for a correct waste disposal.

ABBREVIATIONS

AST : Aspartate Amino Transferase

CLSI : Clinical and Laboratory Standards Institute

CV% : Coefficient of Variation Percentage

EP : Evaluation Protocols
GLP : Good Laboratory Practice

GOT : Glutamate Oxaloacetate Transaminase

IU : International Unit mA : miliabsorbance

mL : mililiter

NCCLS: National Committee for Clinical Laboratory

Standards

QC : Quality Control

REFERENCES

- Young, DS. Effect of Drugs on Clinical Laboratory Tests, AACC Press, 3rd Edition, Washington, (1990).
- Tietz Textbook of Clinical Chemistry, Second Edition, Burtis- Ashwood (1994).
- CCLM 2002; 40(7):725-733, Schumann et al. -IFCC reference procedure for aspartate aminotransferase.
- 4. J. Clin.Chem.Clin.Biochem 8 (1970) 658; 10 (1972) 182
- Clinical and Laboratory Standards Institute (formerly NCCLS). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2004. NCCLS Document EP05-A2.
- International Federation of Clinical Chemistry, Committee on Reference Systems for Enzymes, Chem Clin Lab Med 2002; 40 (7): 725-733.

7. HU Bergmeyer - Methods of enzymatic analysis, (1987).

SYMBOLS

Only for invitro diagnostic use

LOT Lot of manufacturing

R1 Reagent 1

R2 Reagent 2

CONC Concentration

INGRED Reagent Ingredients

REF Reference Number (Catalog No)

SN Serial Number

Expiration date

Storage temperature interval

Read the directions

Biological risk

 ϵ

Archem Diagnostics Industry LTD. ŞTİ.

Organize Sanayi Bölgesi, Mutsan Sanayi Sitesi M8 Blok No: 48 Başakşehir / ISTANBUL TURKEY Tlf: + 90 212 444 08 92

Fax: +90 212 629 98 89 info@archem.com.tr www.archem.com.tr

Rev: V2.3 Date: 01.13 GOT (AST, GOT) Page 3 / 3