IMMUNOGLOBULIN G

(IgG Turbidimetrik)

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KISALTMALAR

CV% : Coefficient of Variation Percentage

GLP : Good Laboratory Practice

IgG : Immunoglobulin G

IU : International Unit

mA : miliabsorbans

mL : mililitre

NCCLS: National Committee for Clinical Laboratory Standards

QC : Quality Control

REFERANSLAR

1. Tietz NW. Clinical Guide to Laboratory Tests. 3rd ed. Philadelphia, PA: WB Saunders Company; 1995:358-359.
2. Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
3. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.
4. Gitlin D, Edelhock H. A study of the reaction between human serum and its homologous equine antibody through the medium of light scattering. J Immunol. 1951;66:67-77.
5. Johnson P, Ottewill RH. A light scattering study of diphtheria toxin-antitoxin interaction. Disc Faraday Soc. 1954;18:327.
6. Marrack JR, Grant RA. The interaction of antigen and antibody in low concentrations of salt. Brit J Exp Path. 1953;34:263.
7. Marrack JR, Richards CB. Light scattering studies of the formation of aggregates in mixtures of antigen and antibody. J Immunol. 1971;20:1019-1038.
8. Narayanan S. Method-comparison studies on immunoglobulins. Clin Chem 1982; 28: 1528-1531.
9. Price CP, Spencer K and Whicher J. Light-scattering immunoassay of specific proteins: a review. Ann Clin Biochem 1983; 20: 1-14.
10. Hellsing K. The effect of hyaluronate, chondroitin

sulfate and chondroitin sulphateprotein complex on the precipitin reaction. Biochem J. 1969;112:475.

11. Hellsing K. Influence of polymers on the antigen-antibody reaction in a continuous flow system. In: Automated Immunoprecipitin Reaction. Colloquium on AIP. Tarrytown, NY: Technicon Inst. Corp.; 1972:17.
12. Hellsing K. Polysacchande - enhanced precipitin reactions with antigens of various sizes. Biochem J. 1969;114:145.
13. Hellsing K, Laurent TC. The influence of dextran on the precipitin reaction. Acta Chem Scand. 1964;18:1303.
14. Hellsing K. The effect of dextran on the reaction between iodine-125 labeled human serum albumin and gamma G-globulin from rabbit anti-albumin sera. Acta Chem Scand. 1966;20:1251.
15. Clinical and Laboratory Standards Institute (formerly NCCLS). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2004. NCCLS Document EP05-A2.
16. Libby RL. New and rapid quantitation technique for the determination of potency of types I and II antipneumococcal serum. J Immunol. 1938;34:269.
17. Goldberg RJ, Campbell DH. The light scattering properties of an antigen-antibody reaction. J Immunol. 1951;66:70.
18. Dati F et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference range for 14 proteins in serum based on the standardization against the IFCC/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
19. Ritchie RF. A simple, direct and sensitive technique for measurement of specific protein in dilute solution. J Lab Clin Med. 1967;70:512-517.

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**Archem Diagnostics Industry LTD. ŞTİ.**Organize Sanayi Bölgesi, Mutsan Sanayi Sitesi