IMMUNOGLOBULIN A

(IgA) - Turbidimetrik

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KISALTMALAR

CV% : Coefficient of Variation Percentage

GLP : Good Laboratory Practice

IgA : Immunoglobulin A

IU : International Unit

mA : miliabsorbans

mL : mililitre

NCCLS: National Committee for Clinical Laboratory Standards

QC : Quality Control

REFERANSLAR

1. Tietz NW. Clinical Guide to Laboratory Tests. 3rd ed. Philadelphia, PA: WB Saunders Company; 1995:354-355.
2. Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
3. Young DS. Effects of Drugs on Clinical Laboratory Tests. 3rd ed. Washington: AACC Press; 1990.
4. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press,1997.
5. Johnson P, Ottewill RH. A light scattering study of diphtheria toxin-antitoxin interaction. Disc Faraday Soc. 1954;18:327.
6. Libby RL. New and rapid quantitation technique

- for the determination of potency of types I and II antipneumococcal serum. J Immunol. 1938;34:269.
7. Ritchie RF. A simple, direct and sensitive technique for measurement of specific protein in dilute solution. J Lab Clin Med. 1967;70:512-517.
8. Dati F et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference range for 14 proteins in serum based on the standarization against the IFCC/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
9. Gitlin D, Edelhock H. A study of the reaction between human serum and its homologous equine antikor through the medium of light scattering. J Immunol. 1951;66:67-77.
10. Narayanan S. Method-comparison studies on immunoglobulins. Clin Chem 1982; 28: 1528-1531.
11. Goldberg RJ, Campbell DH. The light scattering properties of an antigen-antikor reaction. J Immunol. 1951;66:70.
12. Price CP, Spencer K and Whicher J. Light-scattering immunoassay of specific proteins: a review. Ann Clin Biochem 1983; 20: 1-14.
13. Helsing K. Influence of polymers on the antigen-antikor reaction in a continuous flow system. In: Automated Immunoprecipitin Reaction. Colloquium on AIP. Tarrytown, NY: Technicon Inst. Corp.; 1972:17.
14. Helsing K. The effect of dextran on the reaction between iodine-125 labeled human serum albumin and gamma G-globulin from rabbit anti-albumin sera. Acta Chem Scand. 1966;20:1251.
15. Helsing K. The effect of hyaluronate, chondroitin sulfate and chondroitin sulphate protein complex on the precipitin reaction. Biochem J. 1969;112:475.
16. Helsing K, Laurent TC. The influence of dextran on the precipitin reaction. Acta Chem Scand. 1964;18:1303.
17. Helsing K. Polysacchande - enhanced precipitin reactions with antigens of various sizes. Biochem J. 1969;114:145.
18. Clinical and Laboratory Standards Institute (formerly NCCLS). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2004. NCCLS Document EP05-A2.
19. Marrack JR, Richards CB. Light scattering studies of the formation of aggregates in mixtures of antigen and antikor. J Immunol. 1971;20:1019-1038.
20. Marrack JR, Grant RA. The interaction of antigen and antikor in low concentrations of salt. Brit J Exp Path. 1953;34:263.

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**Archem Diagnostics Industry LTD. ŞTİ.**