

Adjustable Pipette

Features

- Adjustable pipette from 0.5 ul to 1000 ul.
- Ergonomic design, comfortable handfeel, simplicity of operation.
- Easy calibration, easy maintenance.
- Digital volume display.
- Built-in ejection of the tip preventing contamination.
- Precise pipette volume.

Technical specifications

Model	Range(ul)	Step volume(ul)	Volume(ul)	Systematic error%	Random error%
SP-10	0.5-10	0.1	5	±1.2	±0.7
SP-50	10-50	0.5	25	±0.5	±0.3
SP-200	20-200	1	100	±0.8	±0.3
SP-1000	100-1000	5	500	±0.8	±0.3

Other model pipettes on request.

sinnowa reserves the right to change the technical specifications without prior notice.
welcome to contact us for latest information.



SINNOWA Medical Science & Technology Co.,Ltd

Add:Qilin industry park,Nanjing,china

Fax:0086-25-84127199

Email:info@sinnowa.com

Tel:0086-25-84127928

P.C: 211135

WWW.SINNOWA.COM